



















# restor3d

## Symbols Glossary

Symbol	Definition	Reference / Standard
	Batch Code	ISO 15223-1: 2016 5.1.5
	Catalog number	ISO 15223-1: 2016 5.1.6
	Quantity	N/A
	Do not re-use	ISO 15223-1: 2016 5.4.2
	Consult Instructions For Use.	ISO 15223-1: 2016 5.4.3
	Consult Instructions For Use. Available on the restor3d website at <a href="http://www.restor3d.com">www.restor3d.com</a>	ISO 15223-1: 2016 5.4.3
	Temperature limitation	ISO 15223-1: 2016 5.3.7
	Keep dry	ISO 15223-1: 2016 5.3.4
	Keep away from sunlight	ISO 15223-1: 2016 5.3.2
	Date of manufacture	ISO 15223-1: 2016 5.1.3
	Manufacturer	ISO 15223-1: 2016 5.1.1
	Use by date	ISO 15223-1: 2016 5.1.4

	Sterilized using ethylene oxide	ISO 15223-1: 2016 5.2.3
	Sterilized using radiation	ISO 15223-1: 2016 5.2.4
	Non-Sterile	ISO 15223-1: 2016 5.2.7
	Do not resterilize	ISO 15223-1: 2016 5.2.6
	For prescription use only	21 CFR 801.15(c)(1)(i)F 21 CFR 801.109
	Do not use if the package has been damaged or opened	ISO 15223-1: 2016 5.2.8